

<i>D</i>	<i>W</i>	<i>S</i>	<i>Method</i>	<i>Analyte</i>	<i>Plas</i>	<i>Gla</i>	<i>Amb</i>	<i>Size</i>	<i>Preserve</i>	<i>Hold AQ</i>	<i>Hold SW</i>
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	API GRO	Gasoline Range Orga	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	40 ML	none	14 d	
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	API PHC	Total Petroleum Hydr	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	250 ML	none	7 d	14 d
<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	EPA 1010	Ignitability	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	100 G	none	28 d	28 d
<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	EPA 1311	TCLP	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 G	none	14 d	14 d
<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	EPA 1312	SPLP	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 G	none	14 d	14 d
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	EPA 160.4	Volatile suspended so	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	none	7 d	
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	EPA 160.5	Settleable solids	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	1 L	none	48 h	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	EPA 1604:MI Agar 40	E.coli	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	120 ML	10 mg Na-thiosulfate; sterile cup/ baq	30 h	
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	EPA 1604:MI Agar 40	E.coli	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	120 ML	10 mg Na-thiosulfate; sterile cup/ baq	6 h	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	EPA 1604:MI Agar 40	Total Coliform	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	120 ML	10 mg Na-thiosulfate; sterile cup/ baq	30 h	
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	EPA 1604:MI Agar 40	Total Coliform	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	120 ML	10 mg Na-thiosulfate; sterile cup/ baq	6 h	
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	EPA 1664A	Oil&Grease/TPH/HE	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	1 L	2 ml H2SO4	28 d	28 d
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	EPA 1666	Pharmaceutical waste	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	40 ML	none	14 d	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	EPA 180.1	Turbidity	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	none	48 h	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	EPA 200.7	Boron	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	500 ml	2 ml HNO3	28 d	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	EPA 200.7	Metals	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	1 L	2 ml HNO3	6 mo	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	EPA 200.7	Metals, dissolved solu	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	1 L	none	ASAP	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	EPA 200.7	Silica, dissolved	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	500 ml	none	28 d	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	EPA 200.7	Silicon	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	1 L	2 ml HNO3	6 mo	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	EPA 200.8	Metals	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	1 L	2 ml HNO3	6 mo	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	EPA 200.8	Metals, dissolved solu	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	1 L	none	ASAP	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	EPA 200.9	Metals	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	1 L	2 ml HNO3	6 mo	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	EPA 200.9	Metals, dissolved solu	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	1 L	none	ASAP	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	EPA 245.1	Mercury	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	1 L	2 ml HNO3	28 d	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	EPA 300.0	Anions	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	1 L	none	28 d	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	EPA 300.0	Bromate	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	500 ML	none	28 d	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	EPA 300.0	Bromide	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	none	28 d	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	EPA 300.0	Chlorate	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	500 ML	none	28 d	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	EPA 300.0	Chloride	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	none	28 d	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	EPA 300.0	Chlorite	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	500 ML	0.5 ml ethylene diamine	14 d	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	EPA 300.0	Fluoride	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	500 ML	none	28 d	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	EPA 300.0	Nitrite	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	none	48 h	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	EPA 300.0	o-phosphate	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	none	48 h	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	EPA 314.0	Perchlorate	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	500 ML	none	28 d	
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	EPA 350.1	Ammonia	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	1 ml H2SO4	28 d	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	EPA 351.2	TKN	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	1 L	2 ml H2SO4	28 d	28 d
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	EPA 353.2	Nitrate/nitrite	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	1 ml H2SO4 in the lab	28 d	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	EPA 420.2	Phenols, total	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	500 ML	1 ml H2SO4	28 d	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	EPA 504.1	DBCP & EDB	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	40 ML	3 mg Na-thiosulfate	14 d	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	EPA 505	Organohalide Pesticid	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	40 ML	5 mg Na-thiosulfate	14 d	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	EPA 515.3	Herbicides	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	2x500 ML	40 mg Na-thiosulfate	14 d	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	EPA 521	UCMR List-2	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	2 x 500 ML	/L:0.08-0.1g Na-thiosulfate	14 d/<10C	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	EPA 524.2	VOCs & TTHMs	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	40 ML	27 mg ascorbic acid + 2 dr 50%HCl	14 d	

<i>D</i>	<i>W</i>	<i>S</i>	<i>Method</i>	<i>Analyte</i>	<i>Plas</i>	<i>Gla</i>	<i>Amb</i>	<i>Size</i>	<i>Preserve</i>	<i>Hold AQ</i>	<i>Hold SW</i>
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	EPA 525.2	SOCs	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	2x1 L	50 mg Na-sulfite + 50%HCl pH<2	14 d	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	EPA 525.2	UCMR List-2	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	2 x 1 L	/L:0.04-0.05g Na-sulfite + 4ML 6N HCl pH<2	14 d/<10C	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	EPA 527	UCMR List-1	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	2 x 1 L	/L:0.1g ascorbic acid+0.35g EDTA-triNa+9.4g KdiHc	14 d/<10C	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	EPA 529	UCMR List-1	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	2 x 1 L	/L:0.5g/L CuSO4+0.3g Tris+4.7g Tris-HCl	14 d/<10C	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	EPA 531.1	SOCs	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	2x250 ML	7.5 ml monochloroacetic acid buffer+20mg Na-thios	28 d	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	EPA 535	UCMR List-2	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	2 x 500 ML	/L:0.1g Ammonium chloride	14 d/<10C	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	EPA 547	Glyphosate	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	40 ML	3 mg Na-thiosulfate	14 d	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	EPA 548	Endothall	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	2x250 ML	20 mg Na-thiosulfate	7 d	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	EPA 549	Diquat	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	1 L	2 ml H2SO4; entire bottle covered with Al-foil	7 d	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	EPA 552.2	Haloacetic Acids	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	250 ML	25 mg NH4Cl	14 d	
<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	EPA 6010B	Metals	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	20 G	none		6 mo
<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	EPA 6020B	Metals	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	20 G	none		6 mo
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	EPA 608	Pesticides & PCBs	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	1 L	none	7 d	
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	EPA 624	VOCs	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	40 ML	none	7 d	
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	EPA 625	SVOCs	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	1 L	none	7 d	14 d
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	EPA 7470	Hg	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	1 L / 20 G	2 ml HNO3	28 d	6 mo
<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	EPA 7471	Hg	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	20 G	none	ASAP	6 mo
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	EPA 8081A	Pesticides	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	1 L / 100 G	none	7 d	14 d
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	EPA 8082	PCBs	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	1 L / 100 G	none	7 d	14 d
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	EPA 8260B	VOCs	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	40 ML / 20 G	none	14 d	14 d
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	EPA 8270C	SVOCs	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	1 L / 100 G	none	7 d	14 d
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	EPA 900 etc.	Radiochemicals	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	4 L	2 ml HNO3	6 mo	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	EPA 901.1	Gamma Spec Scan	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	1 L	2 ml HNO3	180 d	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	EPA 901.1	Iodine 131	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	1 L	none	8 d	
<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	EPA 9010	Cyanide, total	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	none		14 d
<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	EPA 9013	Cyanide, extractable	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	none		14 d
<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	EPA 9014	Cyanide, amenable	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	none		14 d
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	EPA 9020B	TOX (X2)	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	1 ml H2SO4	28 d	
<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	EPA 9023	TOX	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 G	none	28 d	
<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	EPA 9045C	pH	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	100 G	none		ASAP
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	EPA 905.0	Strontium 90	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	1 L	2 ml HNO3	180 d	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	EPA 906.0	Tritium	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	500 ML	none	180 d	
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	EPA 9095	Paint Filter Lq. Test	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	100 G	none	28 d	28 d
<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	EPA/625/R-92/013 pa	Viable Helminth Ova	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	100 G	none; sterile cup/bag		30 d
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Gamma Scintillation	Isotope ID	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	1 L	2 ml HNO3	none	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	HACH 8021	Chlorine, free residual	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	none; collect in Whirl-Pak Bag	ASAP	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	HACH 8167	Bromine	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	none; collect in Whirl-Pak Bag	ASAP	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	HACH 8167	Chlorine, total residua	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	none; collect in Whirl-Pak Bag	ASAP	
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	LAL GEL CLOT	Endotoxin	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	40 ML	none; collect with headspace	7 d	
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	Part 503; app.G	Salmonella	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	120 ML/G	10 mg Na-thiosulfate; sterile cup/ bag	24 h	
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Sabouraud	Yeast/mold in air	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		none; sterile air filter or plate	ASAP	
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	Sec.7.3S	cyanide, reactive	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	100 G	none	14 d	14 d
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	Sec.7.3S	Sulfide reactive	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	100 G	none	14 d	14 d

D	W	S	Method	Analyte	Plas	Gla	Amb	Size	Preserve	Hold AQ	Hold SW
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	SM2120B	Color	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	none	48 h	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	SM2150B	Odor	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	none	ASAP	
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	SM2310B	Acidity	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	250 ML	none	14 d	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	SM2320B	Alkalinity	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	none	14 d	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	SM2340C	Hardness, total	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	1 ml HNO3	6 mo	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	SM2510B	Conductivity	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	none	28 d	
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	SM2510B	Specific conductance	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	none	28 d	
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	SM2540B	Total Residue	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	none	7 d	
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	SM2540B	TS	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	none	7 d	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	SM2540C	Total Dissolved Solid	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	none	7 d	
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	SM2540D	TSS	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	none	7 d	
<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	SM2540G	TVS	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	none	7 d	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	SM2550B	Temperature	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	100 ML	none	ASAP	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	SM3500Ca-B	Calcium as CaCO3	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	2 ml HNO3	6 mo	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	SM4500CN-CE	Cyanide, free	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	300 mg Ascorbic A. & 6 pellets NaOH	14 d	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	SM4500CN-CE	Cyanide, total	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	300 mg Ascorbic A. & 6 pellets NaOH	14 d	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	SM4500CN-G	Cyanide, amenable	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	300 mg Ascorbic A. & 6 pellets NaOH	14 d	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	SM4500H+B	pH	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	250 ML	none	ASAP	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	SM4500NO2-B	Nitrite	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	none	48 h	
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	SM4500O-G	Dissolved oxygen	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	BOD bottle	none	ASAP	
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	SM4500P-E	o-phosphate as P	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	none	48 h	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	SM4500P-E	Total phosphorus	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	1 ml H2SO4	28 d	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	SM4500S2-F	Sulfide	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	3 full droppers (6ml) of 22% zinc acetate	7 d	7 d
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	SM4500SiO2 C	Silica, dissolved	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	500 ml	none	28 d	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	SM4500SO3-B	Sulfite	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	250 ML	none	ASAP	
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	SM5210B	BOD/CBOD	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	1 L	none	48 h	
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	SM5220C	COD	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	1 ml H2SO4	28 d	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	SM5310C	DOC	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	preserve in lab	24 h	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	SM5310C	TOC	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	1 ml H2SO4	28 d	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	SM5540C	Foaming agents MBA	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	none	48 h	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	SM5550B	Tannin & lignin as phe	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	250 ML	none	30 d	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	SM5910B	UV254	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	500 ML	none	48 h	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	SM6040D	Geosmin & MIB	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	250 ml or 4xV	none	14 d	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	SM7500Rn-B	Radon in water	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	2x40 ML	none; no headspace and no label or writing on vial	4 d	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	SM9215B	Heterotropic bacteria	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	120 ML	10 mg Na-thiosulfate; sterile cup/ bag	24 h	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	SM9215B	Pseudomonas	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	120 ML	10 mg Na-thiosulfate; sterile cup/bag; using P.isol.a	ASAP	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	SM9215B	SPC - pour	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	120 ML	10 mg Na-thiosulfate; sterile cup/ bag	24 h	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	SM9215C	SPC - spread	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	120 ML	10 mg Na-thiosulfate; sterile cup/ bag	24 h	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	SM9215D	SPC - filter	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	120 ML	10 mg Na-thiosulfate; sterile cup/ bag	24 h	
<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	SM9221E	Fecal coliform (EC br	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	120 ML	none; sterile cup/bag	24 h	
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	SM9222B	Total Coliform (M EN	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	120 ML	10 mg Na-thiosulfate; sterile cup/ bag	6 h	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	SM9222B	Total Coliform (M EN	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	120 ML	10 mg Na-thiosulfate; sterile cup/ bag	30 h	
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	SM9222D	Fecal coliform (enum	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	120 ML	10 mg Na-thiosulfate; sterile cup/ bag	6 h	

<i>D</i>	<i>W</i>	<i>S</i>	<i>Method</i>	<i>Analyte</i>	<i>Plas</i>	<i>Gla</i>	<i>Amb</i>	<i>Size</i>	<i>Preserve</i>	<i>Hold AQ</i>	<i>Hold SW</i>
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	SM9222D	Fecal coliform (enum	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	120 ML	10 mg Na-thiosulfate; sterile cup/ bag	30 h	
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	SM9230C	Fecal streptococcus	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	120 ML	10 mg Na-thiosulfate; sterile cup/ bag	6 h	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	SM9240D	Sulfate reducing bact	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	120 ML	none; sterile cup/bag	ASAP	
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	SM9260J	Legionella	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	1 L	100 mg Na-thiosulfate; sterile container	ASAP	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	SM9610B	Yeast/mold pour	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	120 ML	none; sterile cup/bag	ASAP	
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	SM9610C	Yeast/mold spread	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	120 ML	none; sterile cup/bag	ASAP	
<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	SM9610D	Yeast/mold filter	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	120 ML	none; sterile cup/bag		ASAP
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	USGS I 1230-85	Hexavalent Chromiu	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	1 L	none	28 d	